




2237 Avenue André Lasquin
74700 SALLANCHES - France
Phone 04 50 58 02 37
Fax 04 50 93 78 60
Email : sales@anthogyr.com
www.anthogyr.com

anthogyr
A global solution for dental implantology

Photo credit : Anthogyr - All rights reserved - 2012-12

SURGICAL DYNAMOMETRIC WRENCH

NEW !



A surgical asset!


anthogyr
A global solution for dental implantology

The instrument that helps you with ACCURACY:

- Adapting your treatment plan
- Supporting your I.L*. protocol decision-making



Surgical dynamometric wrench
Art. Nb. INCCDC



ERGONOMIC

EASY READING
of the tightening torque

RELIABLE TORQUE INDICATION:
0/30/40/60/80 N.cm

REVERSIBLE FUNCTION
to avoid pointless handlings :
PULL - TURN!

COMPATIBLE CONNECTION

with the entire range of
Anthogyr implant screw
wrenches

EXTREMELY SIMPLE TO DISASSEMBLE AND REASSEMBLE:
EASY CLEANING!



*I.L. = Immediate Loading