

SAFETY - EXCELLENCE - RELIABILITY

MONT BLANC[®] LED

SURGERY HANDPIECES



*Your essential
surgical partner!*

anthogyr

A global solution for dental implantology

The **MONT-BLANC®** Surgery line of handpieces was especially designed for surgical applications ; these instruments will exceed your expectations and ease even delicate procedures.

Incorporating the latest technological advancements and anthogyr's manufacturing expertise, these instruments offer uniquely high quality components and exceptional reliability, in an attractive and ergonomic design.

The quality and unrivalled design of the **MONT BLANC®** handpieces will yield greater efficiency in your practice and will be a trusted partner during surgery.

LED



Implantology and surgery motor **IMPLANTEO®** LED

- **LED Technology lighting** for visual comfort.
- **Large touch-screen** : optimum visibility.
- **Intelligence** : data saving on USB memory key, 10 easily configurable programs...
- **Power and precision**: high torque output (80N-cm).

For additional information, visit our website
www.anthogyr.com



1:3



1:1

Maximum safety

The **MONT BLANC®** surgical contra-angle and handpiece were specifically designed for use with your surgery and implantology motor :

- **Effective and precision external irrigation:**
spray can be adjusted in accordance with the length of the drill.
- **No air is projected from the instrument**, thereby avoiding the risk of necrosis, bacterial transmission, subcutaneous emphysema.
- **Instruments can be disassembled for an efficient cleaning!**
- **Thermo-disinfectable.**

LED Technology*

- **«Day light» quality lighting** provides optimal visibility with exceptionally true color.
- **Larger volume of lighting for more comfortable working space.**
- Power 31 000 Lux.
- Unrivalled LED durability.

Optimum ergonomics

The profiled design of the **MONT BLANC®** range provides **ideal working comfort** and avoids hand contraction resulting in a feeling of hand fatigue.

i.coat® coating

Thanks to its uniform surface limiting bacterial retention zones, the **i.Coat® surface treatment optimises the durability of surfaces**. It provides **an ideal grip with a very comfortable feel in the hand**, even during long procedures in a moist environment.

Reliability and power

- **No vibrating feedback** even during difficult use provides accurate and pleasant work!
- Robust stainless steel components contribute to the durability and **longevity of these instruments.**

Can be dismantled for easy maintenance!

The instruments can be **easily and quickly dismantled without tooling**. Designed for effective decontamination and deep cleaning : solution for ever increasing rigorous sterility standards, especially in a hospital environment!



*Only compatible with the **iMPLANTEO® LED** implantology and surgery motor by anthogyr. These surgical instruments are also available without LED.

APPLICATIONS FOR THE **MONT BLANC**® SURGERY LINE:

- Extraction of wisdom teeth
- Osteotomy
- Apical resection
- Bone regularisation
- Preparation of implant site: piercing of cortical bone

SPECIFIC APPLICATIONS FOR THE **MONT BLANC**® 1:3 SURGICAL CONTRA-ANGLE :

- Hemisection
- Cutting bridges / teeth

SPECIFIC APPLICATIONS FOR THE **MONT BLANC**® 1:1 SURGICAL STRAIGHT HANDPIECE :

- Alveolectomy
- Osteosynthesis

	LED	without light	with light
1:1 Ref. 12400XLED ¹ - Straight Handpiece 1:1	●		
Ref. 12400X - Straight Handpiece 1:1		●	
1:3 Ref. 12200XLED ¹ - Contra-angle 1:3	●		
Ref. 12200X - Contra-angle 1:3		●	
Ref. 12200XL ² - Contra-angle 1:3			●

¹ only compatible with the **IMPLANTEO**® LED implantology and surgery motor by anthogyr.
² contra-angles with light only work on a motor with light.

©Anthogyr -C94C - 2011-09

DISCOVER THE **MONT BLANC**® INSTRUMENTS

A complete range for all your needs: implantology, surgery and general dentistry.

The **MONT BLANC**® range, made in France, is the ideal solution for your dental practice.

For additional information, visit our website

www.anthogyr.com



Photos credits : Anthogyr - Atelier du cyclope - All rights reserved

anthogyr

A global solution for **dental implantology**

3495 Piedmont Rd, Building
11, Suite 710
Atlanta GA 30305 - USA
Office (404) 495 - 5982
Mobile (404) 718 - 9689

2237 Avenue André Lasquin
74700 SALLANCHES - France
Phone. +33 (0)4 50 58 02 37
Fax +33 (0)4 50 93 78 60
Email : sales@anthogyr.com

www.anthogyr.com